

# IMMUNO REAGENTS

## 10-9500 Stop Solution for TMB Substrate



### Benefits

- Ready-to-use
- ELISA-reagent
- Sufficient for 10 96-well plates
- Contains 0.5 M H<sub>2</sub>SO<sub>4</sub>

### Product description

Medicago's Stop solution is especially developed for the EC-Blue Enhanced Substrate (TMB) (Article no. 10-9405) and for use in ELISA procedures. The Stop solution is supplied as a ready-to-use suspension.

For laboratory use only.

### Applications

- Stop solution for TMB substrate
- ELISA

### Directions for use

The product is a ready-to-use suspension. Use lab coat, gloves and goggles and handle with care since it contains acid.

### Shipping and storage

The product is shipped at room temperature and is stable for 2 years when stored at 4-8°C. Store in original container, needs to be kept dark.

### Certifications

Medicago's laboratories and manufacturing site in Uppsala are ISO 9001:2008 certified. Each stage of the manufacturing process is controlled and monitored by stringent quality control procedures to guarantee the highest possible quality and lot-to-lot reproducibility.



### Ordering information

Article no.	Product name	Pack size
10-9500-50ml	Stop Solution for TMB Substrate	50 ml

Specifications	Stop Solution for TMB Substrate (10-9500)
Appearance	Clear liquid solution
Storage	4°C to 8°C
Shelf life	2 years